



the ROHO group
shape fitting technology®



NEW! THE ROHO® HARMONY™ CUSHION

LIGHTWEIGHT.

STABLE.

REDUCES DISCOMFORT.

FUNDING FRIENDLY.

The HARMONY™ Cushion from The ROHO Group provides your clients with the lightweight, stable solution they deserve. By combining a lightweight and comfortable ROHO air flotation cushion with a pre-contoured foam base, encased within a moisture resistant cover, your client can maximize their abilities, independence and freedom without the worry of adding unnecessary weight to their mobility system.

Finding a cushion that reduces discomfort, matches your client's shape and addresses today's challenging funding climate, shouldn't be difficult. Think HARMONY from The ROHO Group.

HARMONY CUSHION

LIGHTWEIGHT

Don't limit client function with a heavy wheelchair cushion. Consumers will spend thousands of dollars to shave a few pounds off of a wheelchair, and then disregard how the cushion impacts the weight of the entire mobility system. The client's function can be limited by a heavy mobility system and increases the potential of long-term shoulder damage. At a lean 1.5 lbs. (0.7 kg), the HARMONY Cushion adds minimal weight to your client's mobility solution.

STABILITY

Stability builds a foundation of function. If a client is working to achieve balance or stability they may have difficulty achieving functional goals. The HARMONY cushion is built on an anatomically shaped, foam platform that enhances a client's seated stability and helps maximize their functional potential.

COMFORT

Many times comfort is dismissed as a luxury, but comfort and function go hand in hand. The HARMONY includes a ROHO air flotation cushion positioned over the foam base that provides a shape fitting interface reducing discomfort through enhanced blood flow and achieving longer sitting times.



FEATURES & BENEFITS

ROHO AIR FLOTATION CUSHION:	Reduced Discomfort
PRE CONTOURED FOAM BASE:	Stability/Ease of Transfer
LIGHTWEIGHT:	Increased Mobility & Function
WASHABLE:	Air Flotation Support Pad may be cleaned with soap and water. Moisture-resistant machine-washable cover.

APPLICATIONS

The HARMONY Cushion is most effective for clients who:

- need a surface that reduces discomfort
- require postural support
- are at moderate to low risk for skin breakdown
- are concerned about weight of mobility system

SPECIFICATIONS

INCLUDES:	HARMONY Cushion Base, Air Flotation Cushion Overlay, cushion cover, operations manual, pump, & patch kit
CONSTRUCTION:	Specially formulated polyurethane die-cut foam base, polyvinylchloride air flotation cushion overlay, flame and moisture resistant cover.
SIZES:	Wide Range, to fit chair widths 14" to 20"/30cm to 51cm
WEIGHT:	Varies by size, approximately 1.5 lbs. / 0.7 Kg
WEIGHT LIMIT:	250 lbs. / 113.4 Kg
WARRANTY:	18- Month Limited Warranty

MORE PRODUCTS FROM THE ROHO GROUP



CONTOUR SELECT™
for maximum positioning,
and skin protection



HIGH PROFILE® QUADRO® SELECT®
for skin protection and
more involved
asymmetrical postures

The ROHO Group, Inc. has a policy of continuous product improvement and reserves the right to amend specifications presented in this brochure.



100 North Florida Avenue • Belleville, IL 62221-5429 • USA • www.therohogroup.com
U.S. & Canada: 1-800-851-3449 • 1-618-277-9173 • fax 1-888-551-3449
Outside the U.S. & Canada • 1-618-277-9150 • fax 1-618-277-6518