




The ROHO® HEAL PAD® Cushion
by ROHO, Inc.

The ROHO HEAL PAD provides an adjustable, therapeutic, counterpressure environment which aids in the healing of heel and elbow ischemic ulcers. The unique design, which offers an array of inflatable, multiple height mini air cells, distributes weight and pressures across the contact surface of the patients problem area helping to maintain blood flow to the wound sight, controlling edema, and facilitating the healing of damaged tissues.

The ROHO HEAL PAD is made with three cell heights and when applied, encompasses the heel and elbow area in a cushion of adjustable mini air cells. The HEAL PAD by ROHO, Inc. can be used as a single patient application or on multiple patients. The ROHO HEAL PAD, a proven therapeutic device for all ages and sizes.

1-800-851-3449 ♦ 618-277-9173 ♦ Fax 618-277-9561 ♦ E-mail: crown@ezl.com



AirLITE™



nexus®

CROWN®
THERAPEUTICS, INC.
The Body's Partner



ADAPTOR®



LOW PROFILE®

Features and Benefits

Multiple Mini Air Cell Height: The HEAL PAD offers maximum therapeutic performance for the treatment and prevention of ischemic heel and elbow ulcers by providing a surface which is contoured to fit your patients problem areas.



Adjustable: Shapes to the individual addressing their specific wound care needs.

Ventilation: Individual air cells allow for air circulation between the cells.

Durable Neoprene Rubber: Flame resistant and resistant to urine, oils and medicines. Easily cleaned with soap or warm water and can be steam autoclaved. Product can be easily repaired extending the overall life of the cushion.

Lightweight: The ROHO HEAL PAD weighs approximately twelve (12) ounces making it easy for the patient to transport and move.

24 Month Limited Warranty

Specifications

The ROHO HEAL PAD comes complete with hand inflation pump, repair kit and operating instructions.

Construction—flame resistant neoprene rubber

Weight—Twelve (12) ounces

Cleaning—Use soap and warm water, or steam autoclave

Applications

To be determined by end user, clinician/therapist

CROWN[®]
THERAPEUTICS, INC.

100 Florida Avenue • Belleville, IL 62221-5430 USA

NOTE: All measurements for cushions are approximate from edge to edge.